

360-0238

300 SERIES DIFFICULT INTUBATION TROLLEY FUNCTIONALITY

FUNCIONALIDAD CARRO INTUBACIÓN DIFÍCIL SERIE 300

Simple static trolleys.

Dimensions 640 x 480 mm.

Adequate height, ergonomics: 980 mm
(work surface - floor).

600 x 400 mm ISO drawers. Heights 50, 100, 200 mm.

1 or 2 ISO drawers 50 mm to organize the vials with
injectable medication.

The colour labels are placed on the syringes to avoid
errors and differentiate the medication.

Carros estáticos simples.


Medidas 640 x 480 mm.

Altura adecuada, ergonomía: 980 mm
(superficie de trabajo - suelo).

Cajones ISO 600 x 400 mm. Alturas 50, 100, 200 mm.

1 o 2 cajones ISO 50 mm para la organización de los
viales con medicación inyectable.

Las etiquetas de colores se van colocando en las
jeringuillas para evitar errores y diferenciar la
medicación.



3708 W

364M6 W.12B.14.26LW(R2).32R.36.44L