

360-0240

300 SERIES ANAESTHESIA TROLLEY FUNCTIONALITY

FUNCIONALIDAD CARRO DE ANESTESIA SERIE 300

Double static trolleys.
Dimensions 900 x 630 mm.
Adequate height, ergonomics: 980 mm
(work surface - floor).
400 x 600 mm ISO drawers. Heights 50, 100, 200 mm.
2 ISO drawers 50 mm to organize the vials with
injectable medication.
The colour labels are placed on the syringes to avoid
errors and differentiate the medication.

Carros estáticos dobles.
Medidas 900 x 630 mm.
Altura adecuada, ergonomía: 980 mm
(superficie de trabajo - suelo).
Cajones ISO 400 x 600 mm. Alturas 50, 100, 200 mm.
2 cajones ISO 50 mm para la organización de los viales
con medicación inyectable.
Las etiquetas de colores se van colocando en las
jeringuillas para evitar errores y diferenciar la
medicación.



SYRINGES

JERINGUILLAS

- MEDICATION LABELS DISPENSER
- GLOVE BOX

- DISPENSADOR ETIQUETAS DE MEDICACIÓN
- CAJA DE GUANTES

PROTECTION ON ALL 4 SIDES

PROTECCIÓN 4 LADOS

DISPOSABLE CONTAINER

CONTENEDOR DESECHABLE

IDENTIFICATION

IDENTIFICACIÓN

ERGONOMIC HANDLE

ASA ERGONÓMICA

MEDICATION DIVIDER SET
(Optional)
CONJUNTO SEPARADORES
MEDICACIÓN (Opcional)



3956 G

ERGONOMIC PEDAL
PEDAL ERGONÓMICO

H5/16
H5/20
H5/24