

## 360-0320

### 100 SERIES EMERGENCY TROLLEY FUNCTIONALITY FUNCIONALIDAD CARRO DE EMERGENCIA SERIE 100

Simple extensive mobility trolleys for use in hospitalization units.  
Dimensions 640 x 480 mm.  
Adequate height, ergonomics: 1005 mm (work surface - floor).  
600 x 400 mm ISO drawers. Heights 50, 100, 200 mm.  
Manual locking for drawers with disposable seals.  
This safety system guarantees that the content of the trolley is intact at the critical time of use.

Carros de alta movilidad simples para uso en unidades de hospitalización.  
Medidas 640 x 480 mm.  
Altura adecuada, ergonomía: 1005 mm (superficie de trabajo - suelo).  
Cajones ISO 600 x 400 mm. Alturas 50, 100, 200 mm.  
Bloqueo manual de cajones con precintos desechables. Este sistema de seguridad garantiza que el contenido del carro se encuentra intacto en el momento crítico de utilización.

SERUM  
SUEROS

CARDIAC  
MESSAGE BOARD  
TABLA DE MASAJE  
CARDÍACO

DEFIBRILLATOR  
2 ROTATION AXLES  
DESFIBRILADOR  
2 EJES DE GIRO

ELEVATED PUSH  
HANDLE  
TIRADOR ELEVADO

PROTECTION  
ON ALL 4 SIDES

PROTECCIÓN 4 LADOS

IDENTIFICATION  
IDENTIFICACIÓN

DISPOSABLE  
CONTAINER

CONTENEDOR  
DESECHABLE

GAS BOTTLE  
Ø 120 / 140  
BOTELLA DE GAS  
Ø 120 / 140

SEAL LOCKING  
BLOQUEO PARA  
PRECINTOS

ERGONOMIC  
HANDLE  
ASA  
ERGONÓMICA

MEDICATION DIVIDER SET  
(Optional)  
CONJUNTO SEPARADORES  
MEDICACIÓN (Opcional)

ERGONOMIC PEDAL  
PEDAL ERGONÓMICO



**1296 R** - 1 medication drawer  
- 1 cajón de medicación

**1208 R** - 2 medication drawers  
- 2 cajones de medicación

**H5/16**  
**H5/20**  
**H5/24**

