

## 360-0325

# 100 SERIES EMERGENCY TROLLEY WITH ELEVATED ACCESSORY FUNCTIONALITY FUNCIONALIDAD CARRO DE EMERGENCIA CON ACCESORIO ELEVADO SERIE 100

Simple extensive mobility trolleys for use in hospitalization units.

Dimensions 640 x 480 mm.

Adequate height, ergonomics: 1005 mm (work surface - floor).

600 x 400 mm ISO drawers. Heights 50, 100, 200 mm.

Manual locking for drawers with disposable seals.

This safety system guarantees that the content of the trolley is intact at the critical time of use.

Carros de alta movilidad simples para uso en unidades de hospitalización.

Medidas 640 x 480 mm.

Altura adecuada, ergonomía: 1005 mm (superficie de trabajo - suelo).

Cajones ISO 600 x 400 mm. Alturas 50, 100, 200 mm.

Bloqueo manual de cajones con precintos desechables.

Este sistema de seguridad garantiza que el contenido del carro se encuentra intacto en el momento crítico de utilización.



SERUM

SUEROS

CARDIAC  
MESSAGE BOARD

TABLA DE MASAJE  
CARDÍACO

ELEVATED PUSH HANDLE

TIRADOR ELEVADO

IDENTIFICATION

IDENTIFICACIÓN

GAS BOTTLE  
Ø 120 / 140

BOTELLA DE GAS  
Ø 120 / 140

ERGONOMIC  
HANDLE

ASA  
ERGONÓMICA

ERGONOMIC PEDAL

PEDAL ERGONÓMICO

DEFIBRILLATOR  
2 ROTATION AXLES

DESFIBRILADOR  
2 EJES DE GIRO

PROTECTION  
ON ALL 4 SIDES

PROTECCIÓN 4 LADOS

DISPOSABLE CONTAINER

CONTENEDOR  
DESECHABLE

SEAL LOCKING

BLOQUEO PARA  
PRECINTOS

MEDICATION DIVIDER SET  
(Optional)

CONJUNTO SEPARADORES  
MEDICACIÓN (Opcional)



**1210 R** - 2 medication drawers  
- 2 cajones de medicación

**H5/16**

**H5/20**

**H5/24**