

# MIGHTY 600

## MIGHTY SERIES

ITEM # F-MTY600



U.S. Patent No. 10,792,203



### ABOUT

Meet the Mighty 600, a powerhouse lift that can accommodate almost any user effortlessly. Our strongest lift with a robust 600 pound weight capacity, strength does not compromise versatility for the Mighty 600. With a seat allowing 4" of additional height adjustment, a 360 degree rotation in either direction and the longest reach of any of our lifts, the Mighty provides accessibility solutions with greater water drafts and setbacks up to 44".

The Mighty lifts can also be purchased with pedestal options, giving you even greater wall clearances. Designed for in-ground anchoring, they can also accommodate existing "standard" Scout 2 and Revolution concrete anchors. Complete with rechargeable 24v battery and wall mounted charger, enjoying a day at the pool has never been so easy.

### FEATURES: ⚙️

- ADA Compliant & UL Safety Certified
- 600 lb (272 kg) weight capacity
- 360 degree direct drive rotation
- Wider 22" seat for added comfort
- Adjustable seat height and foot rest
- Stainless steel construction with powder coat finish
- Limited five (5) year warranty
- Includes 24v rechargeable battery and charger
- Ships common carrier, some assembly required
- Made in USA

### ACCENT COLORS 🎨

Available in a variety of accent and powder coat colors



\*Blue seat standard (not a powder coat option)

# strack Liftsysteme

Strack AG  
Ebnatstrasse 125  
8200 Schaffhausen

Tel. 0800 600 500  
Fax 0800 600 900

www.strack.ch  
info@strack.ch

ZSR-Nr. K 7514.14  
NIF-Nr. 366 824  
MWST-Nr. 567 970

Affordable recreation equipment for people of all abilities.

1.888.687.3552 | [aquacreek.com](http://aquacreek.com) | [sales@aquacreek.com](mailto:sales@aquacreek.com)

# MIGHTY 600


## APPLICATIONS

This lift works on in ground and some partially raised pools and spas. To see if this lift works with your pool, visit our website and download the application guide for this product, or contact our friendly sales team. This lift can also be anchored to concrete or paver decks. Choose the appropriate anchor system for your needs. Anchors sold separately.

Standard Anchor



Semi-Recessed Anchor



Paver Anchor



Pedestal



For access to CADs, Revit drawings and more visit our Partner Login page on our website and become an Aqua Creek partner.

## OPTIONS / ACCESSORIES

- Anchoring options for concrete and paver applications
- Cover to protect your investment from the elements\*
- Cycle attachment to convert to exercise bike
- Headrest for comfortable head support
- Solar charging station to conserve energy
- Pedestal options to help clear taller pool and spa walls
- Pull out leg rest for additional wall clearances
- Rescue spine board option
- Transport cart for easy removal from pool side for storage\*\*
- Upper chest strap for added torso stability


\*Color and material options available

\*\*Transport cart not usable with semi recessed anchor

## PREMIUM COVER



## SOLAR CHARGER



Accessorize your lift with a cover to keep it looking like new and a solar charging station for an environmentally friendly power option.



## CUSTOM COLOR OPTIONS

Match your decor with a variety of custom color options.