

PORTABLE PRO POOL 2

PRO POOL SERIES

ITEM # F - P P P 2



U.S. Patent No. D507,769S



ABOUT

No mounts, no mess, no problem. The Portable Pro 2 offers an affordable and compact solution for pool access without having to anchor into your deck. In a single out-and-down motion, you'll find yourself enjoying a day at the pool in seconds with our easy-to-use design that allows users to control their descent with a two-button submersible remote control.

Designed for in-ground pools where the deck is flush with the top of the pool, this lift is easily rolled away from the pool edge when not in use. Along with a rechargeable battery, adjustable foot rest and flip up arm rests, you can choose from a number of add-ons to customize this lift to your specific needs. Looking to customize even more? Call our amazing Customer Service team today to ask about more color options!

FEATURES:

- CE certified
- 350 lb (158 kg) weight capacity
- Includes 24v rechargeable battery and charger
- No anchor required
- Rechargeable battery with wall mount charger
- Stainless steel construction with powder coat finish
- Limited five (5) year warranty
- Ships common carrier, some assembly required
- Made in USA

ACCENT COLORS

Available in a variety of accent and powder coat colors



*Blue seat standard (not a powder coat option)

strack Liftsysteme

Strack AG
Ebnatstrasse 125
8200 Schaffhausen

Tel. 0800 600 500
Fax 0800 600 900

www.strack.ch
info@strack.ch

ZSR-Nr. K 7514.14
NIF-Nr. 366 824
MWST-Nr. 567 970

Affordable recreation equipment for people of all abilities.

1.888.687.3552 | aquacreek.com | sales@aquacreek.com

PORTABLE PRO POOL 2

APPLICATIONS

This lift works with in-ground pools and spas, and does not clear high curbs or walls. To see if it will work for you, visit our website and download the application guide for this product, or contact our friendly sales team.



For access to CADs, Revit drawings and more visit our Partner Login page on our website and become an Aqua Creek partner.

OPTIONS / ACCESSORIES

- Available with either concrete weight plates or sand ballast**
- Cover to protect your investment from the elements*
- Cycle attachment to convert your lift into an aquatic bike
- Headrest for comfortable head support
- Pull out leg rest for additional wall clearances
- Upper chest strap for added torso stability

** Sand is not supplied with the sand ballast version. Customer supplies their own sand.


* Color and material options available

PREMIUM COVER



Accessorize your lift with a cover to keep it looking like new and a solar charging station for an environmentally friendly power option.

SOLAR CHARGER



CUSTOM COLOR OPTIONS

Match your decor with a variety of custom color options.